

ACNS1721 and ACNS1723 Contouring Atlas

Instructions for viewing the two sample cases in MIM

1. Download the two zip files.
2. Extract the files for case 1 and case 2 to different folders.
3. Click on the mimviewer.exe file in each folder to open a MIM viewer with all the files for that case.

Name	Date modified	Type	Size
bin	1/18/2019 7:11 PM	File folder	
MIMviewer.app	1/18/2019 7:11 PM	File folder	
PATIENTS	1/18/2019 7:11 PM	File folder	
autorun.inf	1/18/2019 7:11 PM	Setup Information	1 KB
DICOMDIR	1/18/2019 7:11 PM	File	1,259 KB
mimviewer.exe	1/18/2019 7:11 PM	Application	80 KB
MIMviewer.ini	1/18/2019 7:11 PM	Configuration sett...	1 KB
splash_mimviewer.jpg	1/18/2019 7:11 PM	JPEG image	108 KB